

Growth Factor Measurement in Animals!

Rat- / Mouse-IGF-I ELISA E25
 Rat- / Mouse-IGFBP-2 ELISA E08
 Rat- / Mouse-IGFBP-3 ELISA E031

Mediagnost developed these test systems as tools for the investigation of the IGF-System in mice and rat models e.g. in research and pre-clinical studies



Mus musculus Source: Wikipedia

Mediagnost Rat- / Mouse-IGF-I ELISA (E25)

The IGF-I ELISA (Insulin-like growth factor-I) measures quantitatively the content of murine and rat IGF-I in **serum** and **plasma** sample solutions as well as in **cell culture medium**.

IGF-I is a major biomarker used in diagnosis of growth disorders like idiopathic short stature but also in acromegaly. To elucidate the role of IGF-I in malignancies, short stature or other diseases mice and rats are used as model organisms.

This ELISA allows the reliable and precise measurement of IGF-I in mouse and rat samples.

Assay characteristics E25:

- ✓ highly specific for IGF-I, no crossreactivity to IGF-II
- ✓ single standards: 0.5, 2.5, 6, 12, 18 ng/ml recombinant IGF-I
- ✓ 5% intra-assay variation
- ✓ sensitivity of 0.029 ng/ml
- ✓ linear range from 1:50 - 1:800 dilution
- ✓ low sample volume required: standard dilution of 1:100
- ✓ control sera included

Mediagnost Rat- / Mouse-IGFBP-3 EIA (E031)

The IGFBP-3 ELISA (Insulin-like Growth-Factor Binding Protein-3) is very well suited for determination of IGFBP-3 concentration in **serum** and **plasma** sample solutions of mice and rats.

IGFBP-3 is the one of the key parameters in diagnosis of growth failure. But beside its role in growth physiology IGFBP-3 influences several physiological and pathological processes, e.g. in tumor biology. Therefore several mouse model systems have been established for investigation of IGFBP-3.

This ELISA allows the reliable and precise measurement of IGFBP-3 in mouse and rat samples.

Assay characteristics E031:

- ✓ only small sample volumes necessary (recommended dilution 1:505)
- ✓ low variance: inter- and intra-assay variation $\sigma < 10\%$
- ✓ 7 single standards of recombinant mouse IGFBP-3: 0.39 - 25 ng/ml
- ✓ sensitivity of 0.09 ng/ml
- ✓ linear range from 1:100 - 1:3200 dilution
- ✓ control sera included

Mediagnost MOUSE- / RAT-IGFBP-2 E08

The Mediagnost mouse IGFBP-2 assay provides an excellent tool for quantitative determination of murine and rat IGFBP-2 in **serum** and **plasma** samples. Specially designed for murine research smallest amounts of serum (**2.5 μ l for duplicates**) allow exact determination of IGFBP-2 concentration (recommended sample dilution 1:100).

Assay characteristics E08:

- ✓ 7 single standards (31.25 - 2000 pg/ml) of recombinant mouse IGFBP-2
- ✓ contains antibodies against whole mouse- and rat-IGFBP-2
- ✓ analytical sensitivity of 0.01 ng/ml
- ✓ intra- and interassay variance $< 10\%$
- ✓ control serum included
- ✓ linearity shown for 1:5 – 1:320 dilution

Mediagnost MOUSE- / RAT-IGF-I ELISA E25

Assay Procedure for Double Determinations:

Reconstitution / Dilution of Reagents		
Standards A-E	Reconstitution in Sample Buffer PP	1 ml
Control Serum KS1	Reconstitution in Sample Buffer PP	500 µl
Control Serum KS2	Reconstitution in Sample Buffer PP	500 µl
Wash Buffer WP	dilute in A. dest. (e.g. total volume of 50 ml in a graduated flask and fill up to 1000 ml)	1:20
Sample + Control Sera KS1 and KS2: dilute 1:100 in Sample Buffer PP, mix immediately, incubate at least for 15 min, max. 2h. Use 50 µl for each well in the assay.		
Before conducting the assay equilibrate all reagents to room temperature.		

Assay Procedure for Double Determinations:

Pipette	Reagent	Position
50 µl	Antibody Conjugate AK	in all wells used
50 µl	Sample Buffer PP (blank)	A1 and A2
50 µl	Standard A (0.5 ng/ml)	B1 and B2
50 µl	Standard B (2.5 ng/ml)	C1 and C2
50 µl	Standard C (6 ng/ml)	D1 and D2
50 µl	Standard D (12 ng/ml)	E1 and E2
50 µl	Standard E (18 ng/ml)	F1 and F2
50 µl	Control Serum KS1	G1 and G2
50 µl	Control Serum KS2	H1 and H2
50 µl	Samples	following wells
Cover the wells with the sealing tape.		
Incubation: 1 h at RT, 350 rpm		
5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Enzyme Conjugate EK	each well
Incubation: 30 min at RT, 350 rpm		
5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Substrate S	each well
Incubation: 30min in the dark RT		
100 µl	Stop Solution SL	each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		

Mediagnost Rat- /Mouse-IGFBP-2 ELISA E08

Assay Procedure in Double Determination:

Preparation of reagents:		Reconstitution:	Dilution:
Before assay procedure bring all reagents to room temperature 20-25°C .			
A-G	Standards	in 1 mL Dilution Buffer VP	-
KS	Control Serum	in 250 µL Dilution Buffer VP	1:100 with Dilution Buffer VP
AK	Antibody Conjugate	-	1:100 with Dilution Buffer VP
EK	Enzyme Conjugate	-	1:100 with Dilution Buffer VP
WP	Washing Buffer	-	1:20 with Aqua dest.
Sample dilution: with Dilution Buffer VP 1:100			
Pipette	Reagenzien		Position
100 µL	Dilution Buffer VP (Blank)		A1/A2
100 µL	Standard A (31.25 pg/mL)		B1/B2
100 µL	Standard B (62.5 pg/mL)		C1/C2
100 µL	Standard C (125 pg/mL)		D1/D2
100 µL	Standard D (250 pg/mL)		E1/E2
100 µL	Standard E (500 pg/mL)		F1/F2
100 µL	Standard F (1000 pg/mL)		G1/G2
100 µL	Standard G (2000 pg/mL)		H1/H2
100 µL	Control Serum	(1:100 diluted)	A3/A4
100 µL	Sample	(1:100 diluted)	in the rest of the wells according the
Cover the wells with the sealing tape.			
Sample Incubation: 1 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Antibody Conjugate AK (1:100 diluted)		In each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Enzyme Conjugate EK (1:100 diluted)		In each well
Cover the wells with the sealing tape.			
Incubation: 0.5 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Substrate Solution S		In each well
Substrat S Incubation: 0.5 h in the Dark at 20-25°C			
100 µL	Stopping Solution SL		In jede Vertiefung
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			