



CE

RANDOX

Urinalysis Controls

Third party controls with ranges for a wide variety of test strip manufacturers



Urinalysis Controls

Randox urinalysis controls are specifically designed for use in the quality control of urine test strips and are suitable for use with both manual and automated methods of dipstick analysis. Ranges are provided for an extensive menu of 13 analytes covering both physical and chemical examination of urine specimens.

- 100% human urine
- Two levels available
- Liquid ready-to-use; suitable for Point-of Care Testing
- Open vial stability of 30 days (or 20 immersions) at +2 - 8°C
- True third party control with ranges available for a wide range of test strip manufacturers

New
Now also available in 25ml ready to use dropper bottles



Analytes		
Albumin	hCG	Protein
Bilirubin	Ketones	Specific Gravity
Blood	Leukocytes	Urobilinogen
Creatinine	Nitrite	
Glucose	pH	

Description	Size	Cat. No.
Urinalysis Control Level 1	12 x 12ml	UC5033
Urinalysis Control Level 2	12 x 12ml	UC5034

Our new urinalysis controls offer the same great benefits as before but are now even easier to use. The convenient dropper bottle allows you to simply drop the control material onto the dipstick for instant results and maximum use.

Description	Size	Cat. No.
Urinalysis Control Level 1	6 x 25ml	UC5078
Urinalysis Control Level 2	6 x 25ml	UC5079

RIQAS – External Quality Assessment (EQA)

RIQAS is the world's largest EQA scheme with over 20,000 participants in more than 100 countries worldwide. The RIQAS bi-monthly urinalysis programme covers 14 parameters in total and is suitable for monitoring the performance of both manual and automated methods of urine dipstick analysis.

RANDOX

Randox Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom
 T +44 (0) 28 9442 2413 F +44 (0) 28 9445 2912 E marketing@randox.com I www.randox.com



LT340 MAR 12